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Prismalysat 1,75cal

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Prismalysat 1,75cal

Medical device IIb

Dialysis solution for continual replacement of the renal function (CRRT - Continuous Renal Replacement Therapy) is used by patients suffering from acute renal failure or intoxication with haemodialysable substances.

1. Intended use

Prismalysat 1,75cal is used as dialysis solution in continuous haemodialysis or continuous haemodiafiltration.

2. Caution

Strictly follow the instructions for use according to the Package Leaflet of the product.

The inorganic phosphate concentration should be measured regularly. Inorganic phosphate must be substituted in cases of hypophosphataemia.

During the entire treatment, haemodynamic condition, liquid balance as well as the electrolyte and acid-base balance must be duly monitored.

Only for use with devices intended for CRRT.

3. Warning

Do not use for direct intravenous infusion.

Do not use as substitution solution.

This product is designated for a single use only. It contains a sterile solution. Repeated use of the product may lead to an infection.

Do not use this product after expiry date printed on the box label and on the back side of the bag.

Use only when the overwrap and the bag containing the solution are undamaged.

All seals must be intact. If leakage is discovered, discard the solution immediately since sterility can no longer be assured.

Use only when the solution is clear and free of any particles.

Provided that the solution needs to be warmed up to the body temperature (+37°C), this procedure must be carefully monitored and it is necessary to check that the solution is clear and free of any particles.

Any unused solution must be disposed of.

Reconstituted solution must be used immediately. Other in-use storage times and condition prior to use are the responsibility of the user and should not be longer than 24 hours including the duration of the treatment.

Store at the temperature from +4°C to +30°C.

4. Description

Prismalysat 1,75cal is a double-chamber bag (4750ml and 250ml) made of PP or SA with a total volume of 5000ml. The final reconstituted solution is obtained after opening the peel seal and mixing the contents of both chamber just before use.

Solutions in both chambers are sterile and free of any bacterial endotoxins.

5. Composition

Composition of small chamber before reconstitution:

Calcium chloride dihydrate	5,15	g/l
Magnesium chloride hexahydrate	2,03	g/l
Lactic acid	5,40	g/l
(as (S)-lactic acid solution 90 % w/w)	6,01	g/l
Glucose anhydrous	22,0	g/l
Water for injection ad	1000	ml

Composition of large chamber before reconstitution:

Sodium chloride	6,46	g/l
Sodium hydrogencarbonate	3,09	g/l
Potassium chloride	0,157	g/l
Water for injection ad	1000	ml

Ions and glucose concentration in the dialysis fluid for immediate use are:

Na ⁺	140	mmol/l
K ⁺	2,00	mmol/l
Ca ²⁺	1,75	mmol/l
Mg ²⁺	0,5	mmol/l
Cl ⁻	111,5	mmol/l
HCO ₃ ⁻	32	mmol/l
Lactate (C ₃ H ₅ O ₂)	3	mmol/l
Glucose anhydrous	6,1	mmol/l
Theoretical osmolarity	297	mOsm/l

6. Indications

Dialysis solution Prismalysat 1,75cal is used to treat patients suffering from acute renal failure or intoxication with dialysable substances. The product is indicated for patients with a tendency to hyperkalaemia (high potassium level in the blood).

7. Contraindications

Hypokalaemia

8. Adverse reaction

Low blood pressure (hypotension), especially when too much fluid is removed from the body. Nausea, vomiting, muscle cramps.

9. Interaction with medicinal products and other forms of interaction

The blood concentration of dialysable drugs may be reduced during the treatment. The use of solution may result in alteration of the patient's plasma electrolyte levels.

10. Instructions for use

After mixing up contents of both chambers the bag with the dialysis solution is connected to the dialysate circuit according to instructions for devices designated for CRRT. It is necessary to avoid contaminations of the dialysis solution and any parts being in contact with the dialysis solution. The bag with Prismalysat 1,75cal is disconnected from the dialysate circuit must be disposed of.

Do not use before mixing up both the chambers.


Prismalysat 1,75cal must not be given intravenously or infused into extracorporeal circuit.

Before its use Prismalysat 1,75cal shall be warmed up closely to the patient's body temperature for it not to be substantially decreased patient's temperature.


11. Manipulation

1) Pulling out of the cardboard box


The double-chamber bag containing the dialyzing solution must be pulled out of the cardboard box immediately before use

A) Tear off the perforation by the arrows on both sides of the cardboard box

B) Push through the perforation on the side on both sides of the cardboard box and open the longer side of the cardboard box.

C) Pull out the packaged bag from the and cardboard box (through the longer side of the cardboard box).

2) Opening of secondary packaging




Do not use sharp objects and take care not to damage the primary packaging.

WHEN TAKING THE BAG OUT OF THE SECONDARY PACKAGING, DO NOT GRIP THE TUBES WITH CONNECTORS.




3) Instructions for mixing both chambers

The contents of both chambers must be mixed immediately before use of the dialyzing solution. Before mixing both chambers, check the composition, production batch number and the use by date. Check that the solution is clear. Check the bag and safety outlet stopper are not damaged. The solution must not be used if the packaging or the safety outlet stopper is damaged.

Bag Type A)

Bag Type B)


A) Tip back the small chamber. Press the small chamber so that the solution of the small chamber is concentrated in the middle. Push towards the centre of the small chamber until the peel seal rips.

B) By pressing on the large chamber rip the peel seal completely.


C) Mix the contents of the bag thoroughly.

WHEN THE CONTENTS OF BOTH CHAMBERS HAVE BEEN MIXED, CHECK THAT THE SEAM IS COMPLETELY OPEN, THE SOLUTION IS CLEAR AND COLOURLESS AND NO FLUID IS LEAKING FROM THE BAG.

4) Opening of stopper

5) Breaking the plug in the stopper
After connecting to the tubes of the device, the plug in the stopper will break.

Sample handling with two-chamber bags can be found at http://www.bio-medica.eu/roztoky_navody/

Prismalysat 1,75cal must be used within 24 hours after mixing. Any unused product residue must be disposed of.

12. Expiry date and way of storage

The expiry date of the Prismalysat 1,75cal is indicated on the bag.

The product must not be used after this date. EXP 24 months.

The solution of Prismalysat 1,75cal must be protected from light and stored under the temperatures from +4°C to +30°C.

13. Information for package disposal

We recommend to dispose the package as potentially hazardous waste.

Symbols on the package:

Attention:
Read instruction for use

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